

Qure Medical Drug Master Filed (DMF) Materials List

Compound #	Tomkins Compound	Polymer type	Hardness	Color	Sp. Gr.	Ash %	ΔpH	0.01N I ₂	Application
2002	BO883	Butyl	44±5	Grey	1.17±0.05	31.6-35.6	-0.11	+0.03	Pharmaceutical/Medical
2003	BO175	Butyl	40±5	Red	1.44±0.05	45.2-49.2	+2.60	+1.00	Pharmaceutical/Medical
2010	BO893	CIIR/EPDM Blend	55±5	Grey	1.32±0.05	41.5-45.2	+0.25	+0.01	Pharmaceutical/Medical
2011	BO857	Butyl	42±5	Grey	1.18±0.05	32.2-36.2	+2.20	+0.11	Pharmaceutical/Medical
2012	BO857F2 Strip	Butyl	42±5	Grey	1.18±0.05	32.0-36.0	+1.50	+0.25	Pharmaceutical/Medical
2013	BO848	Butyl	40±5	Grey	1.13±0.05	25.2-29.2	+2.52	+0.61	Pharmaceutical/Medical
2014	BO848H2	Butyl	40±5	Grey	1.15±0.05	25.2-29.2	+2.52	+0.61	Pharmaceutical/Medical
2016	BO880	Halobutyl/EPDM Blend	44±5	Grey	1.20±0.05	32.5-36.5	+0.75	+0.10	Pharmaceutical/Medical
2017	BO860	Butyl/Chloroprene Blend	52±5	Grey	1.41±0.05	50.2-54.2	+1.85	+0.25	Pharmaceutical/Medical
2019	BO878	Halobutyl/EPDM Blend	45±5	Grey	1.15±0.05	28.3-32.3	+0.75	+1.18	Pharmaceutical/Medical
2023	HO852B3	NBR/CIIR Blend	53±5	Grey	1.18±0.05	25.3-29.3	+1.00	+0.20	Pharmaceutical/Medical
2026	BO854	Butyl/Chloroprene Blend	41±5	Grey	1.28±0.05	38.6-42.6	+0.30	+0.10	Pharmaceutical/Medical
2028	BO849	Butyl	40±5	Grey	1.12±0.05	24.2-28.2	--	--	Pharmaceutical/Medical
2030	BO884	Butyl/EPDM Blend	44±5	Grey	1.17±0.05	30.3-34.3	+1.00	+0.10	Pharmaceutical/Medical
2031	BO861	Butyl/Chloroprene Blend	45±5	Grey	1.58±0.05	53.6-57.6	+0.20	+0.10	Pharmaceutical/Medical
2032	O831	Natural Rubber (Latex)	39±5	Grey	1.29±0.05	35.1-39.1	0.40	+0.15	Pharmaceutical/Medical
2035	BO176	Butyl/Polychloroprene	43±5	Red	1.21±0.05	33.9-37.9	+0.50	+0.10	Pharmaceutical/Medical
2036	BO339	Butyl/Chlorobutyl	50±5	Blue	1.38±0.05	49.7-53.7	--	--	Pharmaceutical/Medical
2038	BO155	Butyl	45±5	Red	1.45±0.05	47.3-51.3	+2.40	+0.60	Pharmaceutical/Medical
2039	BO139	Butyl	44±5	Red	1.45±0.05	48.0-52.0	+2.65	+1.00	Pharmaceutical/Medical
2041	BO181	Butyl	42±5	Red	1.27±0.05	35.5-39.5	+2.60	+0.80	Pharmaceutical/Medical
2042	BO168	Polychloroprene	55±5	Blue	1.52±0.05	28.1-32.1	--	--	Pharmaceutical/Medical
2043	BO168	Butyl/Polychloroprene	53±5	Red	1.43±0.05	49.6-53.6	--	--	Pharmaceutical/Medical
2044	SO185A2	NR/SBR Blend	56±5	Red	1.34±0.05	???	+1.30	+0.90	Pharmaceutical/Medical
2046	O843	Natural Rubber/Latex	39±5	Grey	1.30±0.05	35.1-39.1	+0.40	+0.15	Pharmaceutical/Medical
2049	NO429	Polychloroprene	45±5	Green	1.45±0.05	20.1-24.4	+2.30	+0.50	Pharmaceutical/Medical
2052	IO11A2	NR/IR/SBR	37±5	Tan	0.93±0.05	0.61-2.11	+2.40	+0.50	Pharmaceutical/Medical
2053	O830	Natural Rubber (Latex)	39±5	Pink	1.24±0.05	30.7-34.7	+0.40	+0.15	Pharmaceutical/Medical
2054	SO323	NR/SBR Blend	55±5	Blue	1.32±0.05	36.2-40.2	--	--	Pharmaceutical/Medical
2058	O158	Natural Rubber/Latex	46±5	Red	1.56±0.05	50.2-54.2	+0.60	+0.00	Pharmaceutical/Medical
2062	BO178	Butyl/CIIR/EPDM	44±5	Red	1.16±0.05	30.1-34.1	-1.30	+0.50	Pharmaceutical/Medical
2065	BO188	Halobutyl/EPDM Blend	47±5	Red	1.21±0.05	34.6-38.6	+0.60	+0.10	Pharmaceutical/Medical
2066	BO874	Halobutyl/EPDM Blend	42±5	Grey	1.14±0.05	27.2-31.2	+0.75	+0.10	Pharmaceutical/Medical
2067	BO153A2	Butyl	40±5	Red	1.10±0.05	25.0-29.0	+2.30	+0.60	Pharmaceutical/Medical
2076	BO896	CIIR/EPDM Blend	55±5	Grey	1.31±0.05	41.5-45.5			Pharmaceutical/Medical
2080	IO14B10	Polyisoprene	30±5	Tan	0.93±0.05	1.02-1.06			Pharmaceutical/Medical
2081	B1100A	CIIR/Polyisoprene	50±5	Blue	1.30±0.05				Pharmaceutical/Medical
2082	B1100A	CIIR/Polyisoprene	45±5	Blue	1.30±0.05				Pharmaceutical/Medical
2083	BO898	Butyl/Halobutyl	55±5	Grey	1.32±0.05	38.0-42.0			Pharmaceutical/Medical
1309305	1309305	Polyisoprene	30±5	Tan	0.92±0.05				Medical device